



SURGICAL GOWN

Excellence in Hygiene Solutions

1. Scrub hands and dry.
2. Open pack carefully.
3. Hold gown at neck, unfold.
4. Insert arms (keep hands in cuffs).
5. Assistant ties gown at neck.
6. Wear gloves (closed method).
7. Tie inner belt.
8. Pass outer belt, turn, tie front.
9. Keep hands above waist.



Reinforced



Wraparound



Standard



Sterile unless opened or damaged. Store in cool & dry place.



Prakhar Medics

Visit Our Website



Plot No. 07, FIE, Patparganj Industrial Area, Patparganj, Delhi-110092



011-40571547 9871596367



Customercare@prakharmedics.com, pmedindia@yahoo.com



www.prakharmedics.com

LC No. MFG/MD/2021/000329, MFG/MD/2022/000254

ISO. Reg. No. IND/MDQMS/2024/N/07, 24DQNF87